







AUTO TYPE	NPC CODE	IMAGE	DEPTH	PENETRATION	TYPE & SIZE	VOLUME OF BLOOD	PACKING	SUGGESTED TEST AREAS	SUGGESTED PATIENT USE
PINK	FTM065			1.8mm	26g	Low	Box of 100	Ward blood Glucose Monitoring	Diabetes Clinics, Paediatrics, Geriatrics, Home Use
BLUE	FTM066			1.8mm	21g	Medium	Box of 100	Ward blood Glucose Monitoring	Diabetes Clinics, most normal adult skins
GREY	FTM067			2.4mm	21g	Medium/High	Box of 100	Haemoglobin, cholesterol, blood typing, coagulation, & other diagnostic tests	Paediatrics and most normal adult skins
GREEN	FTM069			2.8mm	21g	High	Box of 100	Haemoglobin, cholesterol, blood typing, coagulation, & other diagnostic tests	Thicker adult skins
STERIHEEL ORANGE: PREEMIE	FTM254			0.85 mm	1.75mm incision length	High	Box of 50	Special baby care and maternity units for heel puncture	Preemies
STERIHEEL YELLOW: NEWBORN	FTM253			1.0mm	2.5mm incision length	High	Box of 50	Special baby care and maternity units for heel puncture	Newborns